

Web Portal

<https://help.virtualmedicalcoaching.com>

Download the latest version for updates (Educators Only)



Option to change password, update cookie consent, and select simulation skin tone

Direct link to Self Assessments (might not be available for all sites)

Check that you are signed in as the right user. Sign out whenever you are finished with your session

VIRTUAL MEDICAL COACHING

- ☰ Simulations
- 📊 Reports
- 🛠 Support centre
- ⬇ Download
- ⚙ Preferences
- 🚪 Sign out

Links

- 📄 Self assessments

Welcome **Awet Russom**, not you?
[Click here to sign out.](#)

Simulations

Select a simulation category



Radiography

Assigned simulations

This page varies depending on the simulations assigned to you.

Possible options include Radiography, Birthing, Radiation Safety and Self-Assessment.

Unless self-Assessment is your only product, it will be linked as part of the navigation panel.

Select a simulation and navigate through pre-selectable scenarios, until you reach the simulation launch options.

Use the breadcrumbs at the top to navigate back through the menu.

Reports

Reports

Your reports

Search simulation, session or user

SIMULATION ▾	SESSION ▾	USER ▾	RECORDED ▾	
Radiography	Pelvis - Anteroposterior - Supine	Nicolas Santos	16/08/2023 18:21pm	View report
Radiography	Facial - Lateral - Erect	Xavier Rounsevell	6/07/2023 11:23am	View report
Radiography	Foot - Weight-bearing Dorsoplantar	Xavier Rounsevell	4/07/2023 11:04am	View report

Search to filter

Reports summarise simulations. They can be searched and sorted by simulation, session and user name. Educators have additional access to learner reports.

Reports

Your reports

Learner reports

Search simulation or session

SIMULATION ▾	SESSION ▾	RECORDED ▾	
Radiography	Ankle - Mortise - Supine	31/08/2023 8:35am	View report
Radiography	Sinus - Lateral	25/08/2023 10:52am	View report
Radiography	Cervical Spine - Extension	24/08/2023 9:33am	View report

Sort by type

Learner

Educator